




Self-Guided Demo Tour
Model 2450-EC / 2460-EC Potentiostats
Electrochemistry Lab Systems



Model 2450-EC/2460-EC Potentiostats

Self-Guided Tour


Use this guide to explore what you can do with the Models 2450-EC and 2460-EC Potentiostats. Applications range from Cyclic Voltammetry to Chronoamperometry and graphical viewing of the results.



Topic	What You Will Experience	Page
Introduction to the Models 2450-EC / 2460-EC Potentiostats	Understanding the Models 2450-EC / 2460-EC	3
Exploring the User Interface	Navigating the user interface	10
Cyclic Voltammetry	Using the Model 2450-ECVCUBE to simulate a cyclic voltammetry test	13
Chronopotentiometry	Configuring the Model 2450-EC or 2460-EC for battery charge/discharge using chronopotentiometry	18
Cyclic Voltammetry	Running the cyclic voltammetry test on your electrochemical cell	21
Electrochemistry Tests Available	Discovering what other electrochemistry tests you can run on the Models 2450-EC and 2460-EC	26
Summary	Final summary of the demonstration	27


Models 2450-EC / 2460-EC Potentiostats

The front panel of the user interface




Model 2450-EC Potentiostats

The rear panel of the instrument



Model 2460-EC Potentiostats

The Rear Panel of The Instrument




Models 2450-EC / 2460-EC Potentiostats

How it works

The Source-Measure Unit (SMU)

Stated in the simplest possible terms, a source-measure unit (SMU) instrument integrates the capabilities of a precision power supply (PPS) with those of a high-performance digital multimeter (DMM) in a single instrument. For example, SMU instruments can simultaneously source or sink voltage while measuring current, and source or sink current while measuring voltage. They can be used as stand-alone constant voltage or constant current sources, voltmeters, ammeters, and ohmmeters, and as precision electronic loads. Their high-performance architecture also allows using them as pulse generators, waveform generators, and automated current-voltage (I-V) characterization systems.




Models 2450-EC / 2460-EC Potentiostats

Operation

The Source-Measure Unit (SMU)

SMU instruments simplify capturing the data needed to characterize a wide range of devices with the SMU instruments' built-in DC sweeps, including linear staircase, logarithmic staircase, and custom sweeps.



Models 2450-EC / 2460-EC Potentiostats

Checklist for the tour

Before beginning the tour, please be sure you have the following items:

- Model 2450 or 2460 Instrument
- Power Cord



- Model 2450 or 2460 Quick Start Guide




- Model 2450-ECVCUBE
Battery Demo Fixture

(Included with the
Model 2450-DEMO-KIT)



- Electrochemistry Flash Drive


- Electrochemistry Translator Cable



Models 2450-EC / 2460-EC Potentiostats

Safety Precautions

Before beginning the tour, please be sure you have read the Model 2450 or 2460 Quick Start Guide sections on Safety Precautions before continuing this self-guided tour.



A GREATER MEASURE OF CONFIDENCE

Safety precautions

Observe the following safety precautions before using this product and any associated instrumentation. Although some instruments and accessories would normally be used with nonhazardous voltages, there are situations where hazardous conditions may be present.

This product is intended for use by qualified personnel who recognize shock hazards and are familiar with the safety precautions required to avoid possible injury. Read and follow all installation, operation, and maintenance information carefully before using the product. Refer to the user documentation for complete product specifications.

If the product is used in a manner not specified, the protection provided by the product warranty may be impaired.

The types of product users are:

Responsible body is the individual or group responsible for use and maintenance of equipment, for ensuring that the equipment is operated within its specifications and operating limits, and for ensuring that operators are adequately trained.

Operators use the product for its intended function. They must be trained in electrical safety procedures and proper use of the instrument. They must be protected from electric shock and contact with hazardous live circuits.

Maintenance personnel perform routine procedures on the product to keep it operating properly, for example, setting the line voltage or replacing consumable materials. Maintenance procedures are described in the user documentation. The procedures explicitly state if the operator may perform them. Otherwise, they should be performed only by service personnel.

Service personnel are trained to work on live circuits, perform safe installations, and repair products. Only properly trained service personnel may perform installation and service procedures.

Keithley Instruments products are designed for use with electrical signals that are measurement, control, and data I/O connections, with low transient overvoltage; and must not be directly connected to mains voltage or to voltage sources with high transient overvoltages. Measurement Category II (as referenced in IEC 60664) connections require protection for high transient overvoltages often associated with local AC mains connections. Certain Keithley measuring instruments may be connected to mains. These instruments will be marked as category II or higher.

Unless explicitly allowed in the specifications, operating manual, and instrument labels, do not connect any instrument to mains.

Exercise extreme caution when a shock hazard is present. Lethal voltage may be present on cable connector jacks or test fixtures. The American National Standards Institute (ANSI) states that a shock hazard exists when voltage levels greater than 30 V RMS, 42.4 V peak, or 60 V DC are present. A good safety practice is to expect that hazardous voltage is present in any unknown circuit before measuring.

Operators of this product must be protected from electric shock at all times. The responsible body must ensure that operators are prevented access and/or insulated from every connection point. In some cases, connections must be exposed to potential human contact. Product operators in these circumstances must be trained to protect themselves from the risk of electric shock. If the circuit is capable of operating at or above 1000 V, no conductive part of the circuit may be exposed.

Models 2450-EC / 2460-EC Potentiostats

Basic User Interface Navigation

The Models 2450-EC and 2460-EC Potentiostats are all-in-one instruments with superior ease of use and measurement performance. The Models 2450-EC/2460-EC offer intuitive instrument operation that minimizes the learning curve and enables faster set up to test results.


Objective: See how simple it is to navigate the instrument




- Power up the SMU. Allow SMU to boot up.
- Review the Measure and Source Windows.
- Touch the Measure and Source **Range** buttons to view options.
- Touch the **Source** button to view the numeric keypad. Try entering a value.
- Swipe the lower part of the screen to view the various windows. Explore the options.



- Press the **HOME** key to return to the main screen.



- Explore the various menu icons and their functions. Familiarize yourself with the configuration settings available.
- If you open a window that has a blue bar on the right side of the screen, this means the entire window can be swiped up or down to provide more options. Example:



- Press the **HOME** key when finished.



Pressing the navigation control opens menus and submenus and selects an entered value.

- Rotate the navigation control to select the **MEASURE Range** button. When a button or icon is highlighted in yellow, that button or icon is now in focus and is controllable by the navigation control.
- Press the navigation control to select the button. Rotate the knob to change the ranges.
- Press the navigation control to select the range.
- Try using the navigation control on other screens like the MENU screen.
- Press the **HOME** key when finished.

Models 2450-EC / 2460-EC Potentiostats

Basic User Interface Navigation

The Models 2450-EC and 2460-EC Potentiostats are all-in-one instruments with superior ease of use and measurement performance. The Models 2450-EC/2460-EC offer intuitive instrument operation that minimizes the learning curve and enables faster set up to test results.

Objective: See how simple it is to navigate the instrument

4 VIEWS Screens

The screenshot shows the 'VIEWS' screen with two tabs: 'Graph' and 'Data'. The 'Graph' tab is active, displaying a plot of current (mA) versus time (s). The 'Data' tab shows a 'READING TABLE' with columns for Buffer Index, Time, Reading, and Scale. The table contains 10 rows of data points.

Buffer Index	Time	Reading	Scale
1	11/24 21:02:49.901982	+0.0111 mA	-1
2	11/24 21:02:50.391576	+0.0078 mA	-1
3	11/24 21:02:50.973771	+0.0056 mA	-1
4	11/24 21:02:51.555953	+0.0060 mA	-1
5	11/24 21:02:51.931647	+0.0083 mA	-1
6	11/24 21:02:52.315929	+0.0055 mA	-1
7	11/24 21:02:52.791543	+0.0062 mA	-1
8	11/24 21:02:53.271532	+0.0060 mA	-1
9	11/24 21:02:53.751550	+0.0054 mA	-1
10	11/24 21:02:54.231524	+0.0068 mA	-1

- Select the Data Graph and Sheet view to review your results quickly.
- To view data in the Graph view, press the **MENU** key, then select **Graph** Views. The graph will update if the SMU is continuously running or plot sweep data.
- In the Graph view, explore the Data and Scale tabs to understand how to control what is plotted, how to change axis values, and how to adjust the scaling of the graph.
- To view data in the Sheet view, press the **MENU** key, then select **Sheet** Views.
- Touching any of the data points will bring up more details about that particular measurement.
- Press the **HOME** key when finished.

5 EVENT LOG Screens

The screenshot shows the 'EVENT LOG' legend with a table of icons and their descriptions.

Icon	Description
⚠	A white-outlined triangle means that the system event log has not logged any new events since the last time the event log was viewed.
ℹ	A blue circle means that an informational event message has been logged. The message is for information only. This is used to indicate status changes or information that may be helpful to the user. It also includes commands if the Log Command option is on.
⚠	A yellow triangle means that a warning event message has been logged. This message indicates that a change occurred that could affect operation.
⚠	A red triangle means that an error event message has been logged. This may indicate that a command was sent incorrectly.

- Touch the indicator in the upper right corner of the front-panel display to view the most recent error, warning, or info message.
- Touch the upper right corner of the HOME page to see an event indicator as described above and as shown. Touch this indicator to display the most recent event.
- Select the **Event Log** button to view all messages in the event log.
- To view details about an event, touch any event in the list. This information helps you to resolve errors with the SMU.
- Another way to open the event log is through the main MENU screen. Press the **MENU** key, and under the System column, touch the **Event Log** icon.
- Press the **HOME** key when you have finished.

EVENT LOG Screens

The screenshot shows the 'EVENT LOG' screen with a table of events and a detailed view of a warning event.

Time	Code	Description
11/24 21:04	5084	⚠ Continuous readings unavailable with output off in resistal

The detailed view shows a warning event with the following information:

- Warning 5084
- Continuous readings unavailable with output off in resistance, high-Z, or 4 wire
- Event Log button

Models 2450-EC / 2460-EC Potentiostats

Basic User Interface Navigation

The Models 2450-EC and 2460-EC Potentiostats are all-in-one instruments with superior ease of use and measurement performance. The Models 2450-EC/2460-EC offer intuitive instrument operation that minimizes the learning curve and enables faster set up to test results.

Objective: See how simple it is to navigate the instrument



On the front panel, press the **HELP** key. The SMU has a context-sensitive help system to speed up your learning curve. When an icon or touchscreen button has focus (icon or button is highlighted with a yellow box), that feature or function has context-sensitive help. If there is no specific focus when you press the **HELP** key, overview information for the screen you are viewing displays.

- Rotate the navigation knob until one of the buttons or icons has focus. Press the **HELP** key to display the help.
- If the help display has a blue bar on the right side of the screen, swipe up or down to display more help information.
- Explore the various **MENU** selections and different screens to see more context-sensitive help.
- Press the **HOME** key when finished.

Summary

The user interface of the Models 2450-EC and 2460-EC allows you to set up most instrument functions and features and perform source and measure operations. You can use the touch-screen display to set up the instrument and tests. You can also use the navigation control knob to highlight an item, and then press the control knob to select it.

The major benefits:

- Large, bright, and easy to read front panel
- Context-sensitive front-panel help
- Swipe-screen capability for more measurement information
- User-friendly configurations from the menu view
- Built-in graphing capability to get to your results quickly
- Data sheet that allows you to view source and measure values quickly
- Error logging so you can quickly diagnose the configuration


Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the Model 2450-EC or 2460-EC for cyclic voltammetry.

1 Setting Up



Return the SMU to a default state and set up the graph function:

- Press the **MENU** key.
- Press the System **Info/Manage** icon.
- Press the **System Reset** button and select **OK**.
- Press the **MENU** key.
- Press the Views **Graph** icon.
- Press the **Data** tab.
- Press the **Graph Type** button and select **Scatter/IV**.
- Press the **HOME** key.

2 Setting Up



- Connect the Model 2450-ECVCUBE Demo DUT to the four banana jacks on the front panel. Make sure the cut-out corner of the DUT is closest to the navigation control. If your Demo DUT has two cutout corners, make sure both cutout corners are facing up.

3 Selecting Cyclic Voltammetry Preloaded Test Script



To execute the script saved in the internal memory of the instrument:

- Press the **HOME** key to return to the home screen.
- Touch the **No Script** indicator as shown in the image above. A list of available scripts is displayed.
- Select **Cyclic Voltammetry** from the list.
- The script will begin immediately.

Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the Model 2450-EC or 2460-EC for cyclic voltammetry.

4 Selecting Cyclic Voltammetry Test Script on the Flash Drive



Your Model 2450-EC or 2460-EC demo system came with preloaded electrochemistry test scripts. The demo system also came with a flash drive that contains all of the test scripts and documentation for your reference.

To execute the script from the USB flash drive:


- Insert the USB flash drive into the front-panel USB port.
- At the top of the Home screen, touch the active script indicator. A list of available scripts is displayed.
- Select **usb1/CyclicVoltammetry.tsp**. The script begins to execute immediately.

5 Acquire Eoc




- When the script has executed, it prompts you to indicate if you are ready to measure the open circuit potential (Eoc). If the experiment is set up and measurements are ready to be taken, select **Yes** to begin the test. Select **No** to quit the program.
- With the demo cube inserted, the Eoc will be measured and the result displayed.

- Select **OK**.



6 Define the Potential Sweep Parameters



After the open circuit voltage is measured, the parameters of the potential sweep must be defined: Number of vertices, voltage magnitude, and reference voltage.

- Select **Three** vertices.

Models 2450-EC / 2460-EC Potentiostats


Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the Model 2450-EC or 2460-EC for cyclic voltammetry.

7

Define the potential sweep parameters




- Define the voltage vertices using the following values for the three potentials:
 - Potential 1: **1.2V**
 - Potential 2: **1.8V**
 - Potential 3: **1.2V**
- After each selection, select the **vs. Eref** reference type.




8

Define the Scan Rate & Number of Cycles



- After sweep parameters have been defined, you are prompted to enter the scan rate in units of millivolts per second. Select **OK** to use the default of 25.
- At the next **Enter number of cycles** prompt, enter **2** and select **OK** to accept.




9

Define the Maximum Current Range & Sampling Interval Units



- At the next prompt, select **1mA** for the maximum current range.
- At the next prompt, select **pts/cycle** for the sampling interval units.




Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.


Objective: Discover how quickly you can configure the 2450-EC or 2460-EC for cyclic voltammetry.

10 Define the Number of pts/cycle




- At the next prompt, enter 101 and select **OK**.

11 Begin the Experiment




- You will now be prompted to select whether you want to save the data to a USB flash drive when the experiment is finished. For now, select **No**.



- The final prompt asks you if you are ready to begin the experiment. Select **Yes** to start the test.

12 Graphing the Results



- While the test is running, you will see the SMU plot the results in real time as the test progresses. After two cycles, the resulting voltammogram for the demo cube battery should look like the above image.
- When the test is complete, you will be prompted to either:
 - **Exit**
 - **Repeat**
 - Start a **New** test
- Choose what you what you would like to do.



Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the Model 2450-EC or 2460-EC for cyclic voltammetry.

13

- You can pinch and zoom the Models 2450-EC and 2460-EC touchscreens (just like your smart phone or tablet) on the Graph tab so you can analyze the data. Try zooming in on some of the data on the graph. You can also pinch the touchscreen to zoom out.
- You also can apply graphing cursors on the graph. Along the bottom of the graph screen, slide the bar until you reach the cursor bar.
- Select the button in the right side of the bar so that you can select **Vertical**, **Horizontal**, or **Both** cursors. Using the touchscreen. Adjust the position of the cursors and note how the X and Y values update.
- To return the graph to its original scaling quickly, touch the SmartScale icon in the upper right corner of the graph screen, as shown by the arrow in the image.



14

Graphing the Data

READING TABLE

Buffer Index	Time	Reading	Source
1	11/25 18:49:56.106743	+0.00093 mA	+1.200000 V
2	11/25 18:49:56.590346	+0.00074 mA	+1.212000 V
3	11/25 18:49:57.070318	-0.00061 mA	+1.224000 V
4	11/25 18:49:57.550328	+0.00057 mA	+1.236000 V
5	11/25 18:49:58.030311	-0.00060 mA	+1.248000 V
6	11/25 18:49:58.510293	+0.00062 mA	+1.260000 V
7	11/25 18:49:58.990309	-0.00061 mA	+1.272000 V
8	11/25 18:49:59.470294	+0.00059 mA	+1.284000 V
9	11/25 18:49:59.950280	-0.00061 mA	+1.296000 V
10	11/25 18:50:00.430290	+0.00066 mA	+1.308000 V

- To view data in the Sheet view, press the **MENU** key.
- In the Views column, select **Reading Table**.
- Touch any of the data points to bring up more details about that specific measurement.
- Scroll through the table by sliding your finger along the mini graph in the table.
- Press the **HOME** key when finished.

Models 2450-EC / 2460-EC Potentiostats

Battery Charge and Discharge

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the Model 2450-EC or 2460-EC for battery charge or discharge using chronopotentiometry.

1 Setting Up



Return the SMU to a default state and set up the graph function:

- Press the **MENU** key.
- Press the System **Info/Manage** icon.
- Press the **System Reset** button. Hit **OK**.
- Press the **MENU** key.
- Press the Views **Graph** icon.
- Press the **Scale** tab.
- Press the **X-Axis Method** button and select **All**.
- Press the **HOME** key.

2 Setting Up



- Connect the Model 2450-ECVCUBE Demo DUT to the four banana jacks on the front panel. Make sure the cut-out corner of the DUT is closest to the navigation control. If your Demo DUT has two cutout corners, make sure both cutout corners are facing up.

3 Selecting Chronopotentiometry Preloaded Test Script



To execute the script saved in the internal memory of the instrument:

- Press the **HOME** key to return to the home screen.
- Touch the **No Script** indicator as shown in the image above. A list of available scripts is displayed.
- Select **Chronopotentiometry** from the list.
- The script will begin immediately.

Models 2450-EC / 2460-EC Potentiostats


Battery Charge and Discharge

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the 2450-EC or 2460-EC for battery charge or discharge.

4

Configure Step 1 to Charge at 15 μ A



When the test script executes a number of settings will appear.

- Press **Step Settings**.
- Press Step: SET VALUE.
- Press Current:SET. Change the current to 15 μ A and then press OK.
- Press **Time:5s** and change the duration time to 60 seconds then press OK.
- Press the <<Back button.

5

Configure Step 2 to Discharge at 15 μ A.




Add a second step to discharge the battery

- Add a second step to discharge the battery by pressing More Steps...
- Press Step2: SET VAL.
- Press Current:SET. Change the current to -15 μ A and press OK.
- Press Time:5s and change the duration time to 60 seconds then press **OK**.
- Press <<Back twice to return to main menu.

6

Set the Potential Range



Change the Potential Range to the 20V range by pressing the **20V** button.

Models 2450-EC / 2460-EC Potentiostats


Battery Charge and Discharge

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Discover how quickly you can configure the 2450-EC or 2460-EC for battery charge or discharge.

7

Adjust Measurement Parameters



Change the measurement interval time to 500ms

- From the main screen, press Measurement
- Press Interval
- Change the time to 500ms/pt
- Press <<Back
- From main menu, press Continue

8

Begin the Experiment

When finished, save to USB flash drive?

No

Yes

- You will now be prompted to select whether you want to save the data to a USB Flash Drive when the experiment is finished. For now, hit **No**.

Ready to begin experiment?

No

Yes

- The final prompt asks you if you are ready to begin the experiment. Hit **Yes** to start the test.

9

Graphing the Results



- While the test is running, you will see the 2450-EC or 2460-EC plot the results in real time as the test progresses. First you will see the battery charging up for 60 seconds, and then you will see the battery discharge for 60 seconds.
- When the test is complete, you will be prompted to either:
 - View Statistics
 - Repeat
 - Change Settings
 - Exist
- Choose to Exit. When asked to Exit program, select Yes.

Models 2450-EC / 2460-EC Potentiostats


Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Models 2450 or 2460 in one of these modes.

Objective: Run the cyclic voltammetry test script on a real electrochemical cell

1

Setting Up




Return the SMU to a default state and set up the graph function:


- Press the **MENU** key.
- Press the System **Info/Manage** icon.
- Press the **System Reset** button and select **OK**.
- Press the **MENU** key.
- Press the Views **Graph** icon.
- Press the **Data** tab.
- Press the **Graph Type** button and select **Scatter/IV**.
- Press the **HOME** key.

2

Setting Up



- Prepare your electrochemical cell.
- Connect the translation cable to the Model 2450-EC or 2460-EC as shown above. Note the cable connector markings to insert into the proper banana jacks.
- Connect the other end of the cable to your 2, 3, or 4 terminal cell as shown.




Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Run the cyclic voltammetry test script on a real electrochemical cell

3 Selecting Cyclic Voltammetry Pre-Loaded Test Script




Execute the script saved in the internal memory of the instrument:

- Press the **HOME** key to return to the home screen.
- Touch the **No Script** indicator, as shown in the image above. A list of available scripts is displayed.
- Select **CyclicVoltammetry** from the list. The script will begin immediately.

4 Acquire Eoc




- When the script has executed, it prompts you to indicate if you are ready to measure the open circuit potential (Eoc). If the experiment is set up and measurements are ready to be taken, select **Yes** to begin the test. Select **No** to quit the program.
- With your cell connected, the Eoc will be measured and the result displayed.



- Select **OK**.

5 Define the Potential Sweep Parameters



After the open circuit voltage is measured, the parameters of the potential sweep must be defined: Number of vertices, voltage magnitude, and reference voltage.

- Select Two, Three, or Four vertices for your test.

Models 2450-EC / 2460-EC Potentiostats


Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Run the cyclic voltammetry test script on a real electrochemical cell

6

Define the Potential Sweep Parameters



Define 2, 3, or 4 voltage vertices:

- Use values for the potentials you want to use for your cell.
- After each selection, select the reference type that you want to use.



7

Define the Scan Rate & Number Of Cycles



- After sweep parameters have been defined, you are prompted to enter the scan rate in units of millivolts per second. Enter a scan rate for your test.
- At the next **Enter number of cycles** prompt, enter a value and select **OK** to accept.




8

Define the Maximum Current Range & Sampling Interval Units



- At the next prompt, select the maximum current range you want to use for your cell test.
- At the next prompt, select the sampling interval units.




Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.


Objective: Run the cyclic voltammetry test script on a real electrochemical cell

9 Enter a Value for Sampling Interval Units Chosen



- At the next prompt, enter a value for the sampling interval units you chose in Step 8.
- Select **OK**.

10 Begin the Experiment




- You will now be prompted to select whether you want to save the data to a USB Flash Drive when the experiment is finished. For now, select **No**.



- The final prompt asks you if you are ready to begin the experiment. Select **Yes** to start the test.

11 Graphing the Results



- While the test is running, you will see the SMU plot the results in real time as the test progresses. Depending on your cell chemistry, the voltammogram may look like the above image.
- When the test is complete, you will be prompted to either:
 - **Exit**
 - **Repeat**
 - Start a **New** test
- Choose what you what you would like to do.



Models 2450-EC / 2460-EC Potentiostats

Cyclic Voltammetry

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Objective: Run the cyclic voltammetry test script on a real electrochemical cell

13

- You can pinch and zoom the Models 2450-EC and 2460-EC touchscreens (just like your smart phone or tablet) on the Graph tab so you can analyze the data. Try zooming in on some of the data on the graph. You can also pinch the touchscreen to zoom out.
- You also can apply graphing cursors on the graph. Slide the bottom display bar at the bottom of the graph screen until you reach the Cursor bar, then select the button on the right side of the bar to select **Vertical**, **Horizontal**, or **Both** cursors.
- Using the touchscreen. Adjust the position of the cursors and note how the X and Y values update.
- To quickly return the graph to its original scaling, touch the SmartScale icon in the upper right corner of the graph screen.

Exploring the Graph



14

Graphing the Data

The figure shows a screenshot of the 'READING TABLE' interface. A yellow arrow points to a data point in the table. The table has columns for Buffer Index, Time, Reading, and Source.

Buffer Index	Time	Reading	Source
1	11/05 12:48:23.670575	+0.00021 mA	+0.600000 V
2	11/05 12:48:23.794896	-0.00011 mA	+0.588000 V
3	11/05 12:48:23.914909	-0.00010 mA	+0.576000 V
4	11/05 12:48:24.034908	-0.00010 mA	+0.564000 V
5	11/05 12:48:24.154903	-0.00013 mA	+0.552000 V
6	11/05 12:48:24.274905	-0.00014 mA	+0.540000 V
7	11/05 12:48:24.394871	-0.00019 mA	+0.528000 V
8	11/05 12:48:24.514902	-0.00014 mA	+0.516000 V
9	11/05 12:48:24.634899	-0.00016 mA	+0.504000 V
10	11/05 12:48:24.754894	-0.00014 mA	+0.492000 V

- To view data in the Sheet view, press the **MENU** key.
- In the Views column, select **Reading Table**.
- Touch any of the data points to bring up more details about that specific measurement.
- Scroll through the table by sliding your finger along the mini graph in the table.
- Press the **HOME** key when finished.

Models 2450-EC / 2460-EC Potentiostats

Electrochemistry Tests Available

The Models 2450-EC and 2460-EC Potentiostat instruments come with a set of preloaded electrochemistry tests, including cyclic voltammetry, chronoamperometry, chronopotentiometry, and others. This reduces set up time when you want to use the Model 2450 or 2460 in one of these modes.

Exploring Further

In addition to the Cyclic Voltammetry test, the Models 2450-EC and 2460-EC come preloaded with the following tests that you can explore:

Linear Sweep Voltammetry:

Potential is swept at a user-programmable scan rate between two defined points while current is measured.

Open Circuit Potential:

Measures the cell potential difference between two electrodes with high input impedance as a function of time.

Potential Pulse and Square Wave with Current Measure:

The instrument sources potential at programmable peak and base levels while current is recorded at a user-defined position on the pulse peak level.

Current Pulse and Square Wave with Voltage Measure:

The instrument sources current at programmable peak and base levels while potential is recorded at a user-defined position on the pulse peak level.

Chronoamperometry:

The potential is stepped to a programmed value while the resulting current is measured as a function of time.

Chronopotentiometry:

The current is stepped to a programmed value while the resulting potential is measured as a function of time

Models 2450-EC / 2460-EC Potentiostats

Summary

This concludes the Model 2450-EC / 2460-EC Potentiostat instruments demo tour. The Model 2450-EC / 2460-EC takes a giant step for test and measurement into today's world of state-of-the-art smart technology. Its innovative graphical user interface and advanced touchscreen technology allow intuitive use and minimize the learning curve so that you learn faster, work smarter, and invent more easily. All this combined with Keithley precision and accuracy will enable you to Touch, Test, Invent™ with the new favorite go-to instrument in the lab for years to come.

To learn more, visit:

<http://www.tek.com/keithley> for detailed white papers, application notes, videos, and other materials. For further product demonstrations or to request a quote, please contact your local Keithley authorized distributor.





071347800



© 2016 Keithley Instruments